

NO. <b>263</b>	
MODEL <b>All Kohler EFI Equipped Engines</b>	
DATE ISSUED <b>7/15</b>	DATE REVISED <b>4/18</b>

## EFI Diagnostic Software 25 761 23-S – SOFTWARE UPDATE

UPDATE 11 for Version 2.11 Software

4/2018

**IMPORTANT:** If you are using Version 2.11 Software


If you are using Version 2.11 of the Diagnostic System for Kohler and the displayed version number (in the lower right of the main screen) is less than 2.11.0011 there is an update to the software that should be installed to provide additional diagnostics for ECH440, ECH940, ECH980, and ELH775 engines. This is a **no-charge** update for 2.11 software users.



To download Update 11, enter the following web site into your web browser address line:

<http://kohler.diagsys.com/>

Under the header **Update 11 for Version 2.11 Software**, select "**Click Here to Download this update**" link.



Follow the IMPORTANT instructions as listed on the next screen shown below.

**Update 11 for Version 2.11 Software**  
If you are using Version 2.11 of the Diagnostic System for Kohler and the displayed version number (in the lower right of the main screen) is less than 2.11.0011 there is an update to the software that should be installed to provide additional diagnostics for some engines.  
Click Here to Download this update  
Once you have installed this update, the Version Number displayed in the lower right of the screen when you start the software should say "Version 2.11.0011"  
Read More

Home > Downloads

**Kohler Diagnostic System Version 2.11.0011 Update**  
This update to the Diagnostic Software includes several new engine selections.  
**IMPORTANT!!!**  
First, Download and save this update to your Desktop.  
Then, Right-Click on the installer and select "Run as Administrator" to install the software.  
(This is VERY IMPORTANT).  
This update will only install if you have Version 2.11 already installed and registered on the computer.  
You must be logged into your computer with ADMINISTRATOR RIGHTS or this install will not complete correctly.

Attachment	Size
<a href="#">KDS_2_11_0011.exe</a>	19.13 MB

« Kohler Diagnostic System Ver 2.09 Service Pack 4 up [Reset209.exe](#) »

Diagnostic Systems Associates  
© 2014-2016 Diagnostic Systems Associates, Inc.

## Kohler Diagnostic System Version 2.11.0011 Update

This update to the Diagnostic Software includes several new engine selections.

**IMPORTANT!!!**

First, Download and save this update to your Desktop.

Then, Right-Click on the installer and select "Run as Administrator" to install the software.  
(This is VERY IMPORTANT).

This update will only install if you have Version 2.11 already installed and registered on the computer.

You must be logged into your computer with ADMINISTRATOR RIGHTS or this install will not complete correctly.

Attachment	Size
<a href="#">KDS_2_11_0011.exe</a>	19.13 MB

« Kohler Diagnostic System Ver 2.09 Service Pack 4 up [Reset209.exe](#) »

At the screen to run or save the update (shown below), CLICK on the downward triangle and choose **SAVE AS** to save this update to your desktop.



Right click the installer icon on your desktop and select "Run as Administrator" to install the software.